

Piggyback Construction for Asset Tracking - Secondary Blood Bag Labels

High-Performance Piggyback Material, White Polypropylene

*End Use Applications:
Blood and Plasma Bags*

*For more information on
FLEXcon's
pressure-sensitive film
solutions for Asset Tracking,
contact your local Sales
Representative or our
Product Identification
Business Team at (508)
885-8300.*

Product ID #: FLX069165



Product:

**Flexcon®
PharmCal® BioMark 069165
(PharmCal® P 300 White V-224
200 Poly SC-6 V-246 Spec 50K-8)**

Benefits:

- Print-receptive film printable via thermal transfer and conventional print technologies
- Biocompatible adhesive (container side) mitigates risk of migration through packaging
- Engineered for extreme temperature cycling and water submersion
- Passes centrifuge process
- Backed with a 50 lb. bleached kraft release liner ideal for roll-form converting
- 21 CFR 175.105 compliant*
- ISO10993-5, - 10 (container side) compliant*
- REACH compliant
- 6-month change notification

*FDA CFR 175.105: Adhesives may be safely used as components of articles intended for use in packaging, transporting, or holding food in accordance with the following prescribed conditions. ISO 10993 Biological evaluation of medical devices: Part 5:Tests for in vitro cytotoxicity; Part 10:Tests for irritation and skin sensitization.

“Fit-for-use” testing is recommended under actual application conditions.

Flexcon®
PharmCal® BioMark 069165
(PharmCal® P 300 White V-224
200 Poly SC-6 V-246 Spec 50K-8)

PRODUCT DATA	VALUE	TEST METHOD
Physical Properties		
Thickness (Mils[microns])	Film	3.1 (79) +/- 10%
	Adhesive	0.8-0.9 (20-23) +/- 0.1 (3)
	Liner	2.0 (51) +/- 5%
	Adhesive	1.4-1.5 (36-38) +/- 0.1 (3)
	Liner	3.1 (79) +/- 10%

[Please click here to view the complete product data sheet.](#)

Product Performance and Suitability

All of the descriptive information, the typical performance data, and recommendations for the use of Flexcon products shall be used only as a guide and do not reflect the specification or specification range for any particular property of the product. Furnishing such information is merely an attempt to assist you after you have indicated your contemplated use and shall in no event constitute a warranty of any kind by Flexcon. All purchasers of Flexcon products shall be responsible for independently determining the suitability of the material for the purpose for which it is purchased. No distributor, salesman, or representative of Flexcon is authorized to give any warranty, guaranty, or make any representation in addition or contrary to the above.

Last Modified On: 02/16/2024

